




Ceiling Track Daily Inspections Pre use of the lift

Pre Check List

1. Make sure you are using the Right sling and that the sling is in good working order such as no tears or fraying on both the body of the sling and the loops. There should also be no stretched out or worn down parts again of the body of the sling and the loops
2. Make sure the lift battery is charged prior to using the lift. (*This is indicated by the flashing green light by the Emergency lower down button, as seen below.*)
3. Make sure the Lift is not *Bent/ Crack/ or visibly Damaged* in anyway, do a visual, physical inspection of the lift.
4. Make sure the Bar Hook does **Not** slide in and out without locking in place, and that each of the 3 settings on both sides of the lift is in working order. The bar hooks needs to be at same settings on both side of the lift when carrying out the transfer.
5. Make sure the Gray Belt tension Guard/ Red Belt tension Guard comes out of the lift freely.
6. Make sure the Gray Belt tension Guard/ Red Belt tension Guard are not fraying.
7. Make sure the Gray Belt is Securely in the distance Strap before lifting.
8. Double check sling loops are securely in the hooks on the lift before completely raising.
9. Make Sure the up and down buttons are working on both the Remote and the on the lift.
10. Make sure it is plugged in after use. Make sure the remote is pushed all the way down in the charger, and the wheel locks are on the cart. (Make sure the Emergency power off button is **NOT** pushed in.)



- Each Lift should be cleaned before and after every use. (Use Accel Intervention Wipes or RTU liquid for disinfecting. Keep wet for 1 minute and allow to air dry.)
- If the lift is **not working correctly**, please lock out tag out the lift so it cannot be used. Make a Maintenance Rec, and inform all staff who need to know.