

ARCOR RAIL in USE



A



B

The Arcor Rail is attached to the frame of the bed. The Rail can be positioned parallel to the bed (**A**) to be used by the client to help turn over and move in bed. The Rail also can be used like the armrest of a chair for transfers (**B**): raise Rail up about an inch, turn away from bed, and slot back down into securing mechanism. [Note, to prevent tripping during transfers, the small front wheels of wheelchair needs to point under the chair – moved from position shown in A to that shown in B. For safety, wheel locks are applied for all transfers.]

FLOOR to CEILING POLES / SUPER POLES



These poles are secured between the ceiling and the floor and can be used in any room in the house, i.e. bathroom, bedroom, kitchen or living room. They allow clients to pull to stand and/or to steady self when moving.