

FLOOR LIFT SLINGS



General Purpose Reusable Sling, Padded, Head Support - NSB -100



Ready Car Reusable Sling, Padded, Head Support - NSB-107



Quick Fit Reusable Sling, Padded, Head Support - NSB-600



Single Patient Use, General Purpose Sling - NSB-151





Hygiene Sling, Padded w/Velcro - NSB-300



Walking Sling, Padded - NSB-510

Back Belt Sling w/Buckle - NSB-500



Back Belt Sling w/Velcro - NSB-500



Cardiac Walking Sling - NSB-510-C



Full Body Sling - NSB-200





General Purpose Reuseable Sling, Padded, Head Support - NSB-100



Positioning Sling, Standard, Reusable - NSB-400-S



Positioning Sling, Single Patient Use - NSB-451-U



Barton Transfer Sheet, Standard, Reusable - B-1410



Barton Transfer Sheet, Single Patient Use - B-1450



Positioning Sling, Standard, Reusable - NSB-400-S



Positioning Sling, Single Patient Use - NSB-451-U



Back Belt Sling w/Buckle - NSB-500



Hygiene Sling, Padded w/Velcro - NSB-300



GENERAL SLINGS

General slings are multi-purpose strap-and-loop slings with built in flexibility to achieve leg and body positioning. General slings have equal application on ceiling lifts and floor lifts. The non proprietary strap-and-loop attachments are industry standard and suitable for use on any carry bar designed for strap-and-loop slings.

GENERAL PURPOSE - REUSABLE

The most common and easiest sling to use. Maintains patient in upright sitting position from Bed to chair

NSB-100-LPH	Padded w/Head Support (LARGE)
NSB-100-XPB	Padded w/Head Support (X-LARGE)
NSB-109-2XPB	Bariatric w/Head Support (1000 lb.)

FULL BODY - REUSABLE

Can be used with single or double above knee leg amputation

NSB-200-MPH	Padded w/Head Support (MEDIUM)
NSB-200-LPH	Padded w/Head Support (LARGE)
NSB-200-XPB	Padded w/Head Support (X-LARGE)

CAR LIFT

Specially designed sling to remove patients in an easy controlled compact way out of any size automobile

NSB-107-MPH	Padded w/Head Support (MEDIUM)
NSB-107-LPH	Padded w/Head Support (LARGE)
NSB-107-XPB	Padded w/Head Support (X-LARGE)

QUICK FIT

Designed with 3 inserts in the sling for upper body and head support. Ability to lift patients from the floor and into a wheelchair or bed

NSB-600-MP	Padded (MEDIUM)
NSB-600-LP	Padded (LARGE)

GENERAL PURPOSE - SINGLE PATIENT USE

NSB-151-MPH	Padded w/Head Support (MEDIUM)
NSB-151-LPH	Padded w/Head Support (LARGE)
NSB-151-XPB	Padded w/Head Support (X-LARGE)

EXTENSION LOOPS

Extension loops thread through permanent loops to allow secure field modifications to the length of any strap. Extension loops are also an economical alternative to a secondary carry bar on a portable ceiling lift.

SPECIAL SLINGS

General slings are multi-purpose strap-and-loop slings with built in flexibility to achieve leg and body positioning. General slings have equal application on ceiling lifts and floor lifts. The non proprietary strap-and-loop attachments are industry standard and suitable for use on any carry bar designed for strap-and-loop slings.

HYGIENE

The Hygiene sling allows patient to be transferred to the toilet with ease. Allows for easy removal of patient clothing

NSB-300-MB	Padded w/Buckle (MEDIUM)
NSB-300-LB	Padded w/Buckle (LARGE)
NSB-300-XB	Padded w/Buckle (X-LARGE)
NSB-300-MV	Padded w/Velcro (MEDIUM)
NSB-300-LV	Padded w/Velcro (LARGE)
NSB-300-XV	Padded w/Velcro (X-LARGE)

BACK BELT

Used with the Acure Sit-to-Stand device for raising and lowering patients from the bed, chair or toilet

NSB-500-MB	Padded w/Buckle (MEDIUM)
NSB-500-LB	Padded w/Buckle (LARGE)
NSB-500-XB	Padded w/Buckle (X-LARGE)
NSB-500-MV	Padded w/Velcro (MEDIUM)
NSB-500-LV	Padded w/Velcro (LARGE)
NSB-500-XV	Padded w/Velcro (X-LARGE)
NSB-510-SW	Padded, Walking (SMALL)
NSB-510-C	*Cardiac, Padded, Walking (ONE SIZE)

**Cardiac sling was designed to support patients with chest tubes and catheters*

POSITIONING

This sling is used for re-positioning and turning and for supine transfers. In terms of public awareness, caregiver injuries from patient handling are commonly associated with lift-and-transfer activities, but institutional experience indicates that as much as 50 percent of injuries may be caused by re-positioning and turning. Positioning slings are used primarily with ceiling lifts, but the standard sling may also be used with a floor lift for lifting an fallen resident.

POSITIONING SLING

Multi functioned sling allows caregivers to position patient up in bed, side to side turning, upper body lifting for chest x-rays, lateral transfers from bed to stretcher or bed to bed. The Single Patient Use Sling can be left underneath the patient.

NSB-400-S	Standard w/Straps (700 lb.)
NSB-409-S	w/Straps (1000 lb.)
NSB-451-U	Single Patient Use w/Straps (700 lb.)

TRANSFER SHEET

The Acure Transfer Sheets are a silicone imbedded product that will allow for reduced friction during lateral transfers and the repositioning of patients. The handles on the sheets allow for safe and easy movement of patients up to 1000 lbs. The Acure Transfer sheets are available in a reusable and single patient use model. Use the transfer sheets with the Acure PTS system for an economical No Lift solution.

B1410-151	Standard (36"x84")
B1420-151	Bariatric (46"x84")
B1450-151	Single Patient Use w/Handles (36"x84")
B1460-151	Single Patient Use, no Handles (36"x84")

STANDARD SLING SIZES



XS, XL AND XXL UP TO 600 LBS
ARE AVAILABLE CUSTOM.

SLING BINDING COLORS

SMALL	PURPLE
MEDIUM	GREY
LARGE	BLUE
EXTRA LARGE	BLACK

SLING LOOP COLORS

SMALL	PURPLE
MEDIUM	GREY
LARGE	BLACK

BINDING AND LOOPS ARE COLOR CODED TO FACILITATE SIZE SELECTION AND LEVEL LOADING.

SLING MATERIALS

FLAME RETARDANT POLYESTER
HIGH - STRENGTH NYLON WEBBING
CLOSED - CELL POLYETHYLENE

BODY STRAPS PADDING

ACURE MEDICAL SLINGS DO NOT CONTAIN LATEX.

SLING-TO-EQUIPMENT COMPATIBILITY

SLING TYPE	FLOOR LIFT	CEILING LIFT	SIT-TO-STAND
GENERAL PURPOSE	■	■	
FULL BODY SUPPORT	■	■	
HYGIENE	■	■	■
POSITIONING		■	
BACK BELT			■



SLING-TO-USER COMPATIBILITY

SLING TYPE	ABLE TO STAND	UPPER BODY	USER ABILITY DEPENDENT	HIGHLY INVOLVED	MOST USERS
GENERAL PURPOSE SLING			■		
FULL BODY SUPPORT SLING		■		■	
HYGIENE SLING					■
POSITIONING SLING					■
BACK BELT		■			



SLING DIMENSIONS

	General Purpose			G.P. w/headsupport			Full body support			Hygiene		
	S	M	L	S	M	L	S	M	L	S	M	L
Width	32½"	38½"	44½"	32½"	38½"	44½"	28"	34"	40"	32"	38"	44"
Commode	10"	18"	21"	10"	18"	21"	6½"	10½"	12½"	15½"	21½"	27½"
Length	42½"	45½"	50"	53½"	58½"	65"	49½"	52½"	56½"	37"	38"	39"
Body	22½"	24"	26½"	33½"	37"	41½"	35½"	37½"	39½"	10"	10"	10"

THE INFORMATION IN THIS BROCHURE IS INTENDED AS A GENERAL GUIDELINE. SLING SELECTION IN EACH INDIVIDUAL INSTANCE SHOULD BE DONE IN CONSULTATION WITH COMPETENT MEDICAL AUTHORITY. REFER TO BARTON MEDICAL'S SLING OWNER MANUALS FOR COMPLETE INSTRUCTIONS ON THE USE AND THE CARE OF SLINGS.